


High-Sensitive MAPPS Analysis for High-Confident Immunogenicity Risk Assessment




Elise Pepermans

 IMMUNESPEC


IMMUNESPEC Advanced immunopeptidomics platform



Identification of MHC presented peptides by affinity purification & MS based identification

- Identification of **neoantigens/tumor antigens** for development of **immunotherapy/ precision medicine** 
- Identification of **pathogen** derived antigens for prophylactic **vaccine** development 
- **Immunogenicity assessment** of biotherapeutics (MAPPS) 
 - High-sensitivity immunopeptidomics analysis: maximize number of identified peptides
 - Minimal sample input
 - Semi-automated platform with high throughput capacity
 - Larger screening panels

IMMUNESPEC Advanced immunopeptidomics platform



1 High-outcome purification


State-of-the-art protocol for optimized solubilization • purification • extraction
Semi-automated workflow • HT

2 Sensitive mass spectrometry


Ultraprecise mass spectrometry analysis using a TimsTOF-SCP

3 Extensive data analysis

Optimized data processing for complete immunogenic reports



MAPPS assay. Risk assessment of your biotherapeutic agent.



All protein therapeutics: potential to elicit unwanted immunogenicity (effect on safety, efficacy, PK, PD)

Health authorities: IND submission – requirement thorough immunogenicity risk assessment

→ Preclinical assays for evaluation and mitigation of immunogenicity risk

MAPPS: Measurement of truly presented T cell epitopes

- Immunopeptidomics analysis of antigen presenting cells loaded *in vitro* with the target biotherapeutic
- Pinpointing all the T cell epitopes of the target biotherapeutic (uptake + lysosomal processing + MHC presentation)
- Overview of putative T-cell immunogenic clusters: immunogenic profile of the biotherapeutic
- Vast majority identified peptides self-peptides: high-sensitivity needed not to miss T cell epitopes
- Correlation immunogenicity risk and # presented peptides & # presented clusters

Case Study 1: High-throughput MAPPS analysis of biotherapeutic loaded moDCs

in collaboration with
ImmunXperts
a Q* Solutions Company


60 samples of moDCs (1 M cells)

10 different donors – 6 conditions:

- 4 biotherapeutics (test articles)
- 1 Positive loading control: BetV1
- 1 unloaded control



- Average 8.5K different MHC peptide identifications/sample
- In total 91.029 different MHC peptides identified




- High numbers of presented MHC peptides are identified using only 1M moDCs per sample
- Per donor: reproducible # MHC peptides

Case Study 1: High-throughput MAPPS analysis of biotherapeutic loaded moDCs

in collaboration with
ImmunXperts
a Q Solutions Company

QC: Length distribution of identified peptides

Collaboration ImmunXperts




- Size distribution of identified peptides is conform for MHC-II presented peptides

Case Study 1: High-throughput MAPPS analysis of biotherapeutic loaded moDCs

in collaboration with
ImmunXperts
a Q Solutions Company

QC: MHC-II presented self-peptides in different samples from the same donor



Heat map HG2A – distribution of identified immune peptides per sample (different conditions – same donor)

- MAPPS analysis of samples from same donor: same presentation pattern (self-peptides): reproducibility

Case Study 1: High-throughput MAPPS analysis of biotherapeutic loaded moDCs

in collaboration with
ImmunXperts
a Q* Solutions Company

Positive loading control: BetV1



Distribution of identified immune peptides from BetV1

- 11 to 78 different BetV1 derived peptides/sample
- Total of 218 different BetV1 derived peptides
- Identified peptides in putative immunogenic clusters.

	D1	D2	D3	D4	D5	D6	D7	D8	D9	D10
#Total different peptides	8980	8986	8846	7864	9097	7976	9465	10546	9047	6455
#Different BetV1 peptides	74	33	27	38	69	31	36	78	21	11
% BetV1 peptides	0,8%	0,4%	0,3%	0,5%	0,8%	0,4%	0,4%	0,7%	0,2%	0,2%

- Numerous overlapping MHC-presented peptides are identified, a crucial factor for pinpointing putative immunogenic clusters with high confidence

Case Study 1: High-throughput MAPPs analysis of biotherapeutic loaded moDCs

in collaboration with
ImmunXperts
a Q* Solutions Company

Positive loading control: BetV1



- Identified peptides in putative immunogenic clusters.
- Different HLA genotype: different clusters

Heat map Betv1 – distribution of identified immune peptides per donor

- ⇒ High numbers of presented MHC peptides are identified using only 1M moDCs per sample
- ⇒ Reproducibility between samples
- ⇒ Putative immunogenic regions are identified by multiple peptides: high-confidence immunogenic profile of test article

Case Study 2: Identification of MHC-II presented peptides from marketed biologics


in collaboration with
ImmunXperts
a Q2 Solutions Company

Explorative study to compare immunogenic profile from different biologics

moDCs (<1 M cells)

10 different donors
5 different biologics

- #Total different MHC peptides: 7437 - 13595




- High numbers of presented MHC peptides are identified using < 1M moDCs per sample
- Per donor: reproducible # MHC peptides

Case Study 2: Identification of MHC-II presented peptides from marketed biologics

in collaboration with
ImmunXperts
a Q* Solutions Company

QC: Length distribution of identified peptides



- Size distribution of identified peptides is conform for MHC-II presented peptides


Case Study 2: Identification of MHC-II presented peptides from marketed biologics

in collaboration with


a Q-Solutions Company

QC: MHC Motif Decon Tool: Attribution identified peptides to HLA-DR, HLA-DP and HLA-DQ


Donor B			
DRB1	11:01:01	15:01:01	74.6%
DPB1	04:01:01	04:01:01	12.9%
DQB1	06:02:01	03:01:01	9.2%



Donor E			
DRB1	03:01:01G	04:01:01G	69.7%
DPB1	04:01:01G	04:01:01G	16.2%
DQB1	03:01:01G	02:01:01G	10.1%



Donor I			
DRB1	11:04:01	07:01:01	70.3%
DPB1	02:01:02	04:01:01	16.1%
DQB1	02:02:01	03:01:01	8.7%



- HLA-DR presented peptides: dominant part of all identified MHC-II peptides
- HLA-DP & HLA-DQ presented peptides: subsidiary part of all identified MHC-II peptides

MHC Motif Decon Tool: Morten Nielsen
 (Kaabinejadian et al, 2022)

Case Study 2: Identification of MHC-II presented peptides from marketed biologics


in collaboration with

a Q* Solutions Company

MHC-II Presented peptides from Infliximab

	Total peptides identified (PSM)	Total different peptides	# Heavy chain peptides	# Light chain peptides	% Infliximab peptides
A	18085	9060	18	0	0.20%
C	17155	8476	20	3	0.27%
D	16971	7512	53	23	1.01%
E	16961	7620	36	18	0.71%
H	28257	13516	60	14	0.55%
I	27225	12690	43	16	0.46%

- **Infliximab: mAb with high prevalence ADAs - relatively high numbers of presented peptides**
- **Clusters: in HC CDR2 & CDR3 but most peptides match to HC & LC constant region**




Case Study 2: Identification of MHC-II presented peptides from marketed biologics

in collaboration with
ImmunXperts
a Q* Solutions Company

MHC-II Presented peptides from Nivolumab

	Total peptides identified (PSM)	Total different peptides	# Heavy chain peptides	# Light chain peptides	% Nivolumab peptides
A	18978	9572	5	0	0.05%
B	20068	9931	0	1	0.01%
C	15910	8041	6	0	0.07%
D	17678	7846	7	10	0.22%
E	16747	7437	2	7	0.12%
F	21004	9912	0	0	0.00%
G	21940	10649	0	0	0.00%
H	28305	13581	1	2	0.02%
I	27190	12694	0	6	0.05%
J	24382	11333	1	0	0.01%



- **Nivolumab: mAb with low prevalence ADAs - relatively low numbers of presented peptides**
- **Clusters: in LC CDR2 but most peptides match to HC & LC constant region**

Case Study 2: Identification of MHC-II presented peptides from marketed biologics

in collaboration with
ImmunXperts
a Q* Solutions Company

⇒ Identification of HLA-DR, -DP and -DQ peptides

⇒ HLA-DR presented peptides: dominant part of all identified MHC peptides

HLA-DP & HLA-DQ presented peptides: subsidiary (but non neglectable) part of all identified MHC-II peptides

⇒ Correlation between # identified peptides and #identified clusters in MAPPS assay and immunogenicity incidence & immunogenicity risk

Case Study 3: Identification of MHC-II presented peptides from B cells, myDCs and moDCs


in collaboration with
ImmunXperts
a Qeiosolutions Company

One donor – comparison MAPPS analysis using moDCs vs myeloid DCs & B cells


B-cells (#800 k): BetV1 & Infliximab

Myeloid DCs (#520k): BetV1 & Infliximab

moDCs (#800k) : BetV1, Infliximab, Trastuzumab, Nivolumab, Atezolizumab, ATR-107, unloaded control



different MHC-II presented peptides identified per sample



Western Blot showing MHC-II expression levels in different cell types (Equivalent of 750 cells loaded per lane)


- High sensitivity allows MAPPS analysis of true APCs (B cells #5k – myeloid DC #12k – moDCs #16k)
- Per sample type: reproducible # identified MHC-II peptides
- #Identified peptides correlate with HLA expression levels

Case Study 3: Identification of MHC-II presented peptides from B cells, myDCs and moDCs

in collaboration with
ImmunXperts
a Qeiosolutions Company

Comparison MAPPS analysis of BetV1 loaded cells: B cells / myeloid DC / moDCs

	Total peptides identified (PSM)	Total different peptides	# BetV1 peptides	% BetV1 peptides
myDC	38606	12491	34	0.27%
B cell	14158	5052	11	0.22%
moDC	34959	15781	88	0.56%




- Same clusters identified in B cells & myDCs found in moDCs from same donor (but higher # peptides)
- True APCs (myDCs & B cells) can be used in high-sensitive MAPPS assay

Case Study 3: Identification of MHC-II presented peptides from B cells, myDCs and moDCs

in collaboration with
ImmunXperts
a Qeios Solutions Company

Comparison MAPPS analysis B cells / Myeloid DCs / moDCs



Overlay identified presented self-peptides at protein level in 3 cell types



Heat map HG2A - distribution of identified immune peptides per sample (different cell types- same donor)

- Differences in presented peptide repertoire at protein level
- For proteins from which relatively high # peptides are presented in the three cell types: presentation patterns are highly similar


Case Study 3: Identification of MHC-II presented peptides from B cells, myDCs and moDCs

in collaboration with

a Q* Solutions Company

Comparison MAPPS analysis of BetV1 loaded cells: B cells / myeloid DC / moDCs

Donor X		
DRB1	04:01:01	15:01:01
DPB1	20:01:01	04:02:01
DQB1	06:02:01	03:02:01



MHC Motif Decon Tool
 Morten Nielsen
 (Kaabinejadian et al, 2022)


- Portion of identified DR / DP /DQ peptides highly similar between different cell types from same donor

Case Study 3: Identification of MHC-II presented peptides from B cells, myDCs and moDCs

in collaboration with
ImmunXperts
a Q*Solutions Company

Comparison MAPPS analysis of different biologics loaded on moDCs

	Total peptides identified (PSM)	Total different peptides	# Heavy Chain peptides	# Light Chain peptides	% Biologic peptides
Infliximab	34904	15797	25	21	0.29%
Trastuzumab	33790	15251	5	4	0.06%
Nivolumab	32250	15290	8	2	0.08%
Atezolizumab	34844	15644	16	9	0.16%
ATR107	34330	15666	20	25	0.29%
Unloaded	34846	15753	0	0	0%



- Correlation between # identified peptides and #identified clusters in MAPPS assay and immunogenicity incidence & immunogenicity risk
- High-sensitivity of MAPPS assay: high-confidence immunogenic profile of test article

Case Study 3: Identification of MHC-II presented peptides from B cells, myDCs and moDCs

in collaboration with
ImmunXperts
a Q* Solutions Company

- ⇒ High sensitivity allows MAPPS analysis of reduced number of moDCs and of true APCs
- ⇒ # Identified peptides correlate with HLA expression levels
- ⇒ Correlation between # identified peptides and #identified clusters in MAPPS assay and immunogenicity incidence & immunogenicity risk
- ⇒ Comparison B cells - myeloid DCs - moDCs:
 - Portion of identified DR / DP /DQ peptides highly similar between different cell types from same donor
 - Presentation patterns for same protein are highly similar
- ⇒ High-sensitivity of MAPPS assay: high-confidence immunogenic profile of biotherapeutic

Conclusion

- Correlation between # identified peptides and #identified clusters in MAPPS assay and immunogenicity incidence & immunogenicity risk
- MAPPS assay: majority identified peptides: self peptides
 - ⇒ high-sensitivity needed for high confident immunogenicity risk assessment
 - ⇒ Maximized # identified peptides: certainty about immunogenic profile of biotherapeutic
- High-sensitive, high-throughput immunopeptidomics platform for ultrasensitive MAPPS assay
 - ⇒ Reduced requirement sample input
 - Use limited # cells
 - Less biotherapeutic required for loading cells
 - Large screening panels possible
 - use moDCs / myDCs / B cells / cyno moDCs
 - ⇒ Pinpoint putative immunogenic clusters with great accuracy
 - ⇒ Full overview including DR / DP / DQ presented peptides
 - ⇒ High throughput – fast analysis

High-confident immunogenic profiling via MAPPS for reliable immunogenicity risk assessment

- ⇒ enhanced compound selection & modulation
- ⇒ higher efficacy
- ⇒ safety

IMMUNESPEC

Acknowledgements



Kurt Boonen
Geert Baggerman
Lauren Thijs
Thomas Van Doninck
Lieselotte Van Antwerpen



Sofie Pattijn
Chloé Ackaert
Aurelie Mazy

ImmuneSpec. Meet the team.



Thomas Van Doninck

Geert Baggerman

Elise Pepermans

Kurt Boonen

Lauren Thijs

Lieselotte Van Antwerpen

Pieter-Paul Strybol

SCIENTIST/PM

CSO & CO-FOUNDER

CEO & CO-FOUNDER

CDO & CO-FOUNDER

LAB TECHNICIAN

LAB TECHNICIAN

BUSINESS DEVELOPMENT