


Micro-Flow Imaging and resonant mass measurements
(Archimedes) evaluated for a quantitative differentiation
of protein particles and silicone oil droplets

Andrea Hawe

26th February 2013

EIP Symposium - Munich

Outstanding solutions
for biopharmaceuticals



Requirements & trends for particle characterization

- authorities require more than data for particles $> 10 \mu\text{m}$ and $> 25 \mu\text{m}$ (light obscuration) and visual inspection:

Quantification: 2 -10 μm (type and amount)

Characterization: 0.1 and 1 μm

- increasing need for particle identification for all size classes:
e.g. differentiation **proteinaceous particles** \leftrightarrow **silicone oil droplets**
is highly relevant for products in DCCs or pre-filled syringes


Methods for silicone oil droplet ↔ protein particle

MFI


- flow microscopy combined with digital image analysis
- size range: about 1 – 70 μm
- differentiation of particles (silicone oil and protein) based on image analysis

RMM

- particle size determination by mass measurement
- size range: 100 nm – 1 μm (Nanosensor)
300 nm – 5 μm (Microsensor)
- differentiation of particles based on density (above or below formulation density)




Resonant mass measurements (*ARCHIMEDES*)



- individual particles are weighed in a mechanically resonating microfluidic channel
→ frequency changes depending on particle mass
- calculation of particle size based on assumed density


Protein particles – silicone oil droplets in product




→ artificially generated protein particles and silicone droplets are representative for particles in marketed products

Set-up of the study


Heat stressed IgG



Silicone oil droplets



Preparation of various spiking ratios




Analysis


Micro-Flow Imaging (MFI)
resonant mass measurements (RMM)

Analysis of individual samples

Silicone oil droplets: more in RMM





Protein particles: more in MFI



Possible reasons for differences between RMM and MFI:


- (i) fragmentation of larger silicone oil droplets within RMM to form smaller ones → higher silicone oil numbers in RMM than in MFI
- (ii) possible blockage of RMM by larger protein particles → underestimation of protein particle numbers
- (iii) different measurement principle

MFI and RMM: separate samples ↔ mixture







- study design was controlled by analyzing separate samples (=only silicone oil or protein particles) compared to the mixture of both
- measured concentration for mixed samples = theoretical sum of silicone oil droplets + protein particles
- individual samples can be used for the theoretical concentration in the mixed samples

Customized filter for MFI data





- Samples used for customize filter development
 - silicone oil droplets
 - heat-stressed IgG
- Parameters included into the filter
 - **intensity mean** = mean intensity value over all pixels within one particle
 - **intensity minimum** = the intensity of the darkest pixel of a particle
 - **intensity standard deviation** = differences between higher and lower intensity values within the same particle
 - **aspect ratio** = shape parameter (1 = spherical; 0"= a needle with an infinite length)
- Cut-offs defined at the mean value of 95% confidence intervals between two populations for sizes 2-9 μm based on polynomial fit
- If all four cut-offs were fulfilled \rightarrow particle marked as silicone oil


Customized filter for MFI data





Customized filter for MFI data



silicone oil droplets



protein particles






— customized data filter

- - - "find similar" filter

→ customized filter combines several particle properties and enables more specific cut-offs

MFI: protein ↔ silicone oil

Mixing ratio of 40:60 (protein particles : silicone oil droplets)




→ reliable differentiation for particles > 2 μm with both approaches

MFI: protein ↔ silicone oil at extreme ratios


Mixing ratio of 95:5

A






Mixing ratio of 15:85

B



→ differentiation by customized filter slightly better for particles > 2 μm

RMM: silicone oil \leftrightarrow protein



→ RMM discrimination is very accurate for all ratios

→ higher STDEV when low particle counts are measured or because of increased coincidence and blockage in case of high protein particle concentration

Comparison of MFI and RMM

- MFI and RMM allow to quantify and distinguish silicone oil droplets and protein particles using different measurement principles
 - **RMM** → suitable for the size range of 500 nm to 2-3 μm
 - **MFI** → image analysis reliable starting at a size of 2-4 μm to $\sim 25 \mu\text{m}$ (find similar, customized filter)
 - discrimination by optical evaluation recommended $> 25 \mu\text{m}$
- combination of both techniques to cover the whole size range


Acknowledgements

- Sarah Zölls (Coriolis)
- Daniel Weinbuch (Coriolis)
- Michael Wiggenghorn (Coriolis)
- **Coriolis Pharma Team**

Coriolis' scientific Advisory Board

- Wolfgang Friess (LMU)
- Gerhard Winter (LMU)
- Wim Jiskoot (Leiden University)
- John Carpenter (University of Colorado)

Outstanding solutions
for biopharmaceuticals



Thanks for your attention



Coriolis Pharma

Biopharmaceutical Research and Development Service

Coriolis Pharma


Biopharmaceutical Research and Development Service

Am Klopferspitz 19 | 82152 Martinsried | Germany

Phone +49 (0) 89 – 5 42 44 98-0 | Fax +49 (0) 89 – 5 42 44 98-22

andrea.hawe@coriolis-pharma.com

Outstanding solutions
for biopharmaceuticals



Cutting edge science for your biopharmaceutical

science-driven formulation
and lyophilization
development tailored for
your product

expert team of more than
20 scientific employees
(pharmacists, engineers,
biochemists, chemists...)

board of expertise
Prof. Dr. Wolfgang Frieß
Prof. Dr. Wim Jiskoot
Prof. Dr. Gerhard Winter
Prof. Dr. John Carpenter



Coriolis Pharma

Biopharmaceutical Research and Development Service

innovative analytical and
technical equipment,
focus aggregate and particle
characterization

cutting edge research in the
field of protein sciences
with top publications

protein, mABs, peptides,
vaccines, oligonucleotides,
colloidal systems, NCEs

www.coriolis-pharma.com

Key expertise and service of Coriolis Pharma

Formulation research and development for biopharmaceuticals

Analytical method development

cGMP analytics for particle characterization

Stability testing (ICH) and stress testing

Lyophilization cycle development, optimization and scale-up

Sterile manufacturing of pre-clinical material